Symplicity® Catheter Instructions For Use [EN] Symplicity® Katheter Bedienungsanleitung [DE] Gebruiksaanwijzing voor de Symplicity®-katheter [NL] Cathéter Symplicity® Mode d'emploi [FR] Catetere Symplicity® Istruzioni per l'uso [IT] Catéter Symplicity® Instrucciones de uso [ES] Cateter Symplicity® Instruções de Utilização [PT] Καθετήρας Symplicity® Οδηγίες Χρήσης [EL] Symplicity®-kateter Brugsanvisning [DA] Symplicity[®] kateter Bruksanvisning [SV]

US: Caution – Investigational Device Limited by United States law to investigational use.

sources of pacing and defibrillation available during ablatio

tor and recommended accessories.

• Do not move the catheter during energy delivery.

Do not advance the Symplicity Catheter while the tip is de-

• The catheter should be deployed while using fluoroscopy.

Avoid treating in arteries with diameter less than 4mm.

• Remove all guidewires from treatment site before activating the

Avoid treating in arteries with significant disease or with flow-

should be provided to the patient. In clinical studies of the

Do not allow the catheter electrode to come into contact with a

metal instrument or surface during energy delivery. Superficial

Radio frequency surgery uses high-frequency output. Do not

perform procedures if flammable or explosive media are pres-

• Safe use of monopolar radio frequency surgery demands prope

surface of the dispersive electrode (also known as neutral

apply the dispersive electrode where fluid may pool.

sheeting is recommended for this purpose.

electrode, return electrode, or grounding pad) may result in a

should be placed on the thigh or other non-bony area out of the

The patient should not come into contact with metal parts which

example operating table supports, etc.). The use of antistatic

are grounded or have an appreciable capacitance to grou

it, such as flammable anesthetics, skin preparation agents, or

coupling of the dispersive electrode to the patient. Follow all of

the manufacturer's directions for the placement of the dispersive

t least 250 seconds was maintained during RF delivery.

and performing a complete analysis of the implanted device

• The catheter is intended for use only with the Symplicity Genera-

LBL039 Rev B © 2010 Ardian, Inc.

Symplicity® Catheter Instructions For Use [EN] PRODUCT OVERVIEW

The Symplicity® Catheter (REF RDN006) is designed to deliver low-level radio frequency (RF) energy through the wall of the renal artery to achieve renal denervation. The Symplicity® Catheter is part of the Symplicity® Catheter System™, which also includes a radio frequency (RF) Generator (the Symplicity® Generator) and foot switch. The catheter is placed in the vasculature adjacent to the target neural site using standard interventional catheter techniques. RF energy is then delivered through the catheter to

The Symplicity Catheter is a 5Fr compatible, single-use RF probe with an effective length of 108cm. It features a monopolar plati electrode at the distal tip of the catheter that is used in conjunction with a standard dispersive electrode (also known as neutral elec trode, return electrode, or grounding pad). The platinum electrode is radiopaque which assists in the positioning of the catheter using scopic guidance. The distal end can be actively deflected using the handle lever. The deflectable section is designed to assist in positioning the electrode within the artery, and the flexibl tip is designed to facilitate consistent wall contact. The catheter is torquable, such that the user can easily rotate the catheter to treat different locations around the vessel. An integrated exte cable is used to connect the catheter to the Symplicity Generator. he catheter is intended for use only with the Symplicity Generator. Figure 1 includes a diagram of the Symplicity Catheter.

The Symplicity Generator delivers radio frequency energy suf

front panel displays information such as temperature, impedance, count down time, and messages. Controls include action utton(s), an RF activation button, and an RF activation foot switc The Symplicity Generator includes an automated algorithm that Catheter. Please refer to the User's Manual for further information on the Symplicity Generator. INTENDED USE The intended use of the device is to deliver low-level radio

frequency (RF) energy through the wall of the renal artery to denervate the human kidney. CONDITIONS FOR USE

 Prior to use, inspect packaging for damage or breach of sterile of damage or breach. • The catheter has a limited shelf life. The product must be used

n or before the "Use By Date" provided on the package. Prior to use, the product should be stored in a cool, dry place. e product should not be exposed to organic solvents, ionizing radiation, or ultraviolet light.

• The Symplicity Catheter is supplied in a sterile, single pouched package in a labeled product carton. • The Symplicity Generator, with foot switch, is non-sterile and

reusable. The Symplicity Generator is supplied separately. Intervention Risks: Death, cardiopulmonary arrest, heart rhythm

or embolism (which may result in ischemic events such as myo-REFERENCES cardial infarction, stroke, kidney damage, or peripheral ischemia physician should consult recent literature on current medical actices for peripheral interventions. damage, arterial spasm, arterial stenosis, arterial dissection or

burns could occur.

bio-intestinal gases.

metal surfaces.

skin burns, and thermal injury to the vasculature or other structure The Symplicity Catheter is designed to work specifically with from energy application. Contrast agents, narcotics, anxiolytics, the Symplicity Generator and should only be used with Ardiar other pain medications and anti-vasospasm agents are also co nufactured or recommended devices. Closely follow the monly utilized during the procedure or after the procedure, and

nstructions for Use and the User's Manual provided with the use of these agents is associated with known risks. Other Treatment Risks: Proteinuria, hematuria, electrolyte dis Attaching the Dispersive Electrode turbances, worsening renal function, hypotension, hypertension repare the patient using standard techniques for electrosurger

nsure the patient's entire body, including extremities, is insulated

Biological Hazards: Risks of infection, toxicity, adverse hematolfrom contact with grounded metal parts. A disposable dispersive ogy, allergy, hemorrhage, and pyrogenicity. olies with IEC 60601-2-2 is required. Closely follow instructions Environmental Hazards: Consistent with standard hospital protocols for proper use and disposal of biological wastes. to achieve good skin contact by the entire adhesive surface of the dispersive electrode may result in a burn or measured high imped-Radiation Hazards: Consistent with normal use of x-ray during ance. The dispersive electrode should be placed on the thigh or interventional procedures. ther non-bony area out of the angiogram field (shave if necessary WARNINGS & PRECAUTIONS r good contact). Do not apply the dispersive electrode where

> Please review the RISKS and WARNINGS & PRECAUTIONS prior Symplicity Generator Setup and Operation 1. If desired, insert the foot switch connector into the receptacle

the generator User's Manual.

eter from the packaging.

the catheter must not be used.

nected to the Symplicity Generator.

movement of the catheter and handle.

nsertion and Positioning of the Catheter

ndicator messages or warnings (e.g., Fault or status lights).

. Additional instructions for using the generator are provided in

Using aseptic technique, carefully remove the Symplicity Cath-

Inspect the catheter for damage. If the catheter is damaged,

3. Pass the cable connector out of the sterile field to be con-

NOTE: The cable should be secured to the table or drape

using a towel clamp, hemostats, or equivalent to help prevent

catheterization is necessary before using this device. • The catheter is intended for single patient use only - Do not resterilize or reuse. Reuse, reprocessing or resterilization may on the front panel. empromise device integrity and may create the risk of trans

Following a system self test, the system is in Standby mode, no • The Symplicity Generator should be powered on and allowed to complete system self-test prior to placing the catheter within only once after the Symplicity Generator powers up. tinue. The user may then press the RF button or depress the foot switch at any time to activate the RF output.

 The Symplicity Catheter System has not been evaluated in patients who are pregnant, nursing, or planning to become

A thorough understanding of the technical principles, clinical ap-

plications, and risks associated with percutaneous transluminal

sion of infectious diseases from one patient to another, which

may result in injury, illness or death of the patient.

• The Symplicity Catheter System has not been evaluated in patients with Type I diabetes mellitus. Avoid use of the Symplicity Catheter System in individuals

Carefully read all instructions prior to use.

in whom a reduction in blood pressure would be considered azardous, e.g., those with hemodynamically significant valvula heart disease. • The Symplicity Catheter System has not been evaluated in

patients with prior renal angioplasty, indwelling renal stents or aortic grafts. tors (ICDs) may be adversely affected by RF ablation. Consider deactivating ICDs during ablation, having temporary external

4. **VERIFY:** That the system components are functional by confirming mit dem Symplicity-Katheter programmiert. Nähere Informationen that the Symplicity Generator displays a valid temperature (e.g., room temperature or approximately 23°C) or message telling user

5. Using standard interventional technique, gain access to the femoral artery and place a 6Fr (or larger) introducer sheath following the manufacturer's instructions. 6. Using standard interventional technique, carefully advance an • Prüfen Sie vor der Verwendung die Verpackung auf Beschädigung appropriate guide catheter into the target renal artery using oder Aufbruch der sterilen Verpackungsversiegelungen. Verwende

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Treat distally in the vessel first and move the catheter proximally for subsequent treatments. Under fluoroscopic guidance, position the treatment electrode at the target treatment site in the vessel.

7. Advance the Symplicity Catheter through the guide catheter to

fluoroscopic guidance.

NOTE: For precise positioning, use the shaft (not the handle) to move the catheter within the vessel. NOTE: If desired, torque the catheter to various orientations. The andle rotator can be used to rotate the shaft separately from the handle body. The handle rotator features tactile clicking every 45 degrees of rotation, and can rotate up to 2 full revolutions in each

NOTE: If desired, deflect the catheter using the handle lever while holding the catheter handle steady. NOTE: If desired, monitor the impedance display to assist in the assessing for sufficient wall contact.

Once intended location is achieved, depress the foot switch or press the RF button to deliver RF energy. The Symplicity will cease power delivery when the treatment is completed. T timer begins counting down and the blue indicator light turns o ntinuously when the foot switch is depressed or the RF button is pressed. At any point in the procedure, the treatment may be stopped by depressing the foot switch or pressing the RF button. $\label{eq:control}$ NOTE: Ensure electrode stability prior to and during treatmen

the foot switch or RF button to stop the Symplicity Generator. The treatment electrode may then be repositioned. NOTE: The Symplicity Generator may automatically stop delivering RF energy if certain conditions are met. A system indicator message or code may appear on the display (see Symplicity or User's Manual) and can be cleared by press CONTINUE button. In case of a hardware fault condition, the

NOTE: If you need to reposition the treatment electrode, depress

light, sound an audio alert, and display a fault code on the display. See Symplicity Generator User's Manual. 4. Once the treatment is completed, reposition the catheter as above to locate the next target treatment site (at least 5mm proximal to

released by pushing the handle lever forward.

NOTE: Ensure that guide catheter is flushed with heparinized Whenever flushing the guide catheter, wait at least 3 seconds before initiating next treatment.

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methods described above and deliver treatments as described

Removing the Catheter 1. Before withdrawing the Symplicity Catheter, ensure the electrode tip is not deflected by pushing the handle lever forward com-

 Increased vessel reactivity, such as spasm, may be encountered 2. Withdraw the catheter completely from the patient. Check the when treating in areas with reduced blood flow, such as dua distal tip of the catheter for coagulum. renal arteries or beyond significant renal artery branch points

3. Withdraw the guide catheter from the patient. • Avoid occluding renal blood flow during energy delivery. 4. Remove the introducer sheath from the femoral artery and use During fluoroscopy, care should be taken to avoid excessive standard of care procedures to achieve hemostasis at the puncture exposure of the patient to contrast agents.

SYMBOLS USED IN LABELING Symplicity Catheter System, an activated clotting time (ACT) of Catalog Number • In clinical studies of the Symplicity Catheter System, individual

LOT Batch Code effort was taken to not re-cross fresh treatment sites. Treatment Use by Date locations spaced < 5mm apart have not been evaluated. STERILE R Contents sterilized by gamma irradiation • In the event that the generator stops due to high impedance or

Attention: See Instructions for Use applied. Use only sterile saline and gauze pad to clean the tip Contents sterile unless package has been damaged or • Ensure that guide catheters are flushed with heparinized saline Do not use if package is damaged.

 Interference produced by the operation of high-frequency surgical equipment may adversely influence the operation of othe Single use only electronic medical equipment such as monitors and imaging Do not resterilize

Hazardous electrical output. This equipment is for use only by EC REP EU Authorized Representative qualified medical personnel trained in the use of this equipment. Manufacture trode at the same time during energy delivery. Superficial skin

DISCLAIMER OF WARRANTY

NOTE: Although the Symplicity Catheter, hereafter referred to as duct." has been manufactured under carefully controlled co tions, Ardian and its affiliate, hereafter referred to as "Ardian," have rol over conditions under which this product is used. Ardi therefore, disclaims all warranties, both expressed and implied with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Ardian shall not be liable to any person or entity for any medical expenses o defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind Ardian to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by court of competent jurisdiction, the validity of the remaining portion of this Disclaimer of Warranty shall not be affected.

PATENTS AND TRADEMARKS

registered trademarks of Ardian, Inc.

The marks ARDIAN, SYMPLICITY as well as the Ardian Logo are

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rade varumärken som tillhör Ardian, Inc
Ardian, Inc. 1380 Shorebird Way Phone: +1 650.417.6534 Fax: +1 650.417.6599
KRAUTH medical KG (GmbH & Co.) Wandsbeker Königstraße 27-29